



Prima Care

HOMECARE & CLINIC FURNITURE





INTRODUCTION

For more and more South Africans, caring for family members at home is the preferred alternative to hospitals or nursing care. With the increasing availability of affordable home care products, it's possible for patients to recover and maintain their health in the privacy of their own home.

However, the broad range of today's home care products can be overwhelming.

To assess your needs and help you choose the product best for you, we have Home HealthCare Consultants with many years of experience in the home health care and nursing fields.

Further we only use high quality products from trusted and well known German and European suppliers.



Our locations are staffed with certified fitters who offer private services in our stores or in the convenience of your own home.

We offer a wide variety of products and services in the following areas:

- + Anti bedsores mattresses
- + Bariatric patients
- + Bathroom safety
- + Care beds
- + Convalescent care
- + Daily living
- + Orthopedic supports
- + Recovery after surgery
- + Walking aids
- + Wheelchairs and accessories

In addition, we understand that home care can be a temporary situation. Therefore, we have rental services for the following products:

- + Electric home care beds
- + Anti bedsores mattresses
- + Overbed tables with locker



ADVANTAGES HOME CARE BED

Electric beds typically use low voltage motors to control and elevate the lying surface of the bed and adjust the height of bed. This is done by a small handset which is used to control independent low voltage motors integrated into the bed.

The choice in electric beds for people with disabilities in homecare and community use has grown considerably over the last year. An electric bed can be adjusted to a number of different positions using a profiling lying surface, which can commonly incline the upper body and raise the lower body independently of each other.

For users with restricted mobility or disabilities which impair movement, the electrically controlled profiling facilitates independent positioning in order to sit up in bed to read or watch TV, which can help to maximise comfort and assist relaxation. Profiling can also be used to sit-up the user prior to transfer into or out of bed. Our adjustable beds offer a chair or comfort position which is perfect for those that want to eat, read or relax in bed.

Our beds have separate head, centre and foot sections which are all independently moveable. The handset controls low voltage electric motors which change the profile of the bed. Powered by a motor underneath the mattress surface, each section of the bed can be positioned independently.

Lowering the height of the bed is beneficial for those with a limited mobility when getting into and out of bed.

The variable height also aids those that transfer from a wheelchair to set the bed at the correct level to assist the manoeuvre into or out of bed. A height adjustable bed can be raised to a high height, which means that carers don't have to bend dangerously when assisting a user. Height adjustability reduces the risk of back injury to the carer and is a requirement by many social services authorities which provide home care.

Some adjustable electric beds can carry patients up to 300 kg. For these specially designed beds, please visit our bariatric section.?



HOMECARE BED FULLY ELECTRIC



Our German home care beds stand for low weight and easy handling. Four different lying surface and side rail combinations are available. In most cases you do not need to worry about heightening pieces for side rails – even if an alternating pressure system is in use. In addition the new SMPS 29V technology is used, meaning that power will only be used when the remote control is pressed. These beds are solid, robust, well-priced and German quality.

TECHNICAL DATA:

Lying surface (W x L)	90 x 200 cm
Dimensions (W x L)	104,5 x 219 cm
Height adjustment	From 40 to 80cm (without mattress)
Safe load	185 kg
Warranty	2 years
Made in	Germany



HEMOCARE BED ACCESSORIES



1



8



2



9



3



10



5

1. Maximum variability by 4-sectional lying surface (metal grid or wooden slats).
2. Height, back rest and foot section electrically adjustable, foot section with knee kink which is adjustable per hand control while lying in bed without nursing – patented!
3. On both sides integrated, twisting proof side rails.
4. Easy assembly and disassembly by one person.
5. Power disconnection: this safety system disconnects the bed from the main power when the hand control is not in operational use.
6. Ergonomically designed hand control with large keys, with integrated lock function by nursing key to prevent misuse.
7. Incl. lifting pole and triangle grip

OPTIONAL EXTRAS:

8. Heightening for side rails, e.g. for use of high cell mattresses
9. Upholstery for side rails
10. Lengthening piece for bed and also for mattress of 20 cm for tall patients



PATIENT LIFTS FULLY ELECTRIC



The original among patient hoists

The Foldy is one of the most sold patient hoists in the home care field. The approved technology is improved constantly. The Foldy now has a removable battery box and thus it is optimized for use in the ambulatory field. The large accessories of slings, scales and spreader bars make the Foldy one of the most flexible devices in this operation area.

You can identify a professional hoist manufacturer by the quality and assortment of the accessories, among others. With us you have a wide range of spreader bars and slings so that full care is ensured perfectly. We offer our slings in different sizes. We assemble your hoist according to the patient's needs.



BARIATRIC CARE

UP TO 300 KG

A bed designed for heavier patients needs a sturdy, reinforced frame, and of course a larger sleeping area. In positioning the bed, it is also important to make sure that there is sufficient room around it to accomplish any nursing tasks or transfers of the patient safely and comfortably. A clearance of 5 feet (1.5 metres) at the sides and foot of the bed allows for positioning of bariatric hoists and wheelchairs and for carers to assist.

For pressure relief and management of respiratory and circulatory problems, often associated with obesity, the profiling frame can be coordinated with a bariatric alternating and lateral rotating mattress, which turns the patient and ensures that pressure is not placed continuously on the same areas, thus freeing the carer from a significant task.



A double scissor construction and 8 castors supports you to solve the problem to care for very heavy patients. This bariatric bed can bear patients up to 300 kg.

Our bariatric range comes in sizes from 90 x 190 cm up to 140 x 240 cm.



HOSPITAL BED ACCESSORIES

European quality hospital beds, fully refurbished and imported from Germany. All beds are suitable for home care as well as hospital use. They allow the care giver to assist the patient with their dignity and comfort as a priority whilst minimizing the effort needed to do so.

The height adjustment allows the patient to easily be removed or placed on the bed with ease, or allow themselves to do the same.

All beds have a movable backrest function that allows it to lift and recline, which allows the user/patient to sit up and lie down with the greatest of ease.

Certain beds have further functions and/or electric to optimize the comfort and movement of the user/patient, promoting blood circulation, independence and one's self-esteem.



Our Extra's are there to make life that much easier: bedside lockers with an over-bed table is a valued combination to finish the bed with a convenient touch. Or even cot sides to ensure the user/patient safety, preventing one won't fall out of bed.



WHAT IS BEDSORE ? WHICH MATTRESS ?

Just about every family in South Africa has a relative that is bedridden for most likely an indefinite period of time. Although this is unfortunate for a number of reasons, the primary concern for many people is preventing bedsores.

Bedsores are lesions that form anywhere on the body due to unrelieved pressure, temperature and several other factors. Fortunately, there are now mattresses specifically manufactured to prevent bedsores from growing.

When the time comes to invest in a proper anti bed sore mattress for your relative, there are 3 important tips that you need to keep in mind:

1. Find a Quality Mattress within Your Budget

Although this seems like a no-brainer, some believe that in order to have a name brand mattress, they have to make great sacrifices with their finances for the sake of their loved one's comfort. Don't fall into this trap. Please keep in mind that there are several variations of anti bed sore mattresses that greatly shift in price, which brings us to our next tip.

2. Know Your Bedsores and Mattresses

There are four distinct types of bedsores you need to be aware of. Stages I-II are the beginning stages of bedsores, while stages III-IV are seen as deep wounds. Today, there are mattresses that are customized to heal bedsores in all stages.

3. Keep the Stage IV Mattress in Mind

The stage IV bed sore mattress is the most expensive model. If you can afford the investment, then it's recommended more than the other models. It's very powerful and has incredible longevity. So, if you want the best for your loved one, the stage IV mattress is what you need.

Bedsores have always been a problem for bedridden individuals. Luckily, there are ways to prevent and treat them. Consider all the types of anti bed sore mattresses mentioned above if you have a relative in this situation – for their assurance and comfort.




PrimaCare has invested towards mattresses of only the best quality, designed for those dependent on their beds, with phenomenal results. Mattresses, designed to reduce pressure on the skin for the greatest protection and comfort are now available in South Africa.

Formulated with advance high memory foam that allows the full support at low pressures for those who spend extended times in their bed.

Many have suffered from soft tissue, nerve and circulation damage, bed sores and as well as pain and discomfort when using ordinary unforgiving mattresses. Resulting in a list of additional expenses and a sluggish recovery from relevant health issues.

We take pride in our superior technology and offering our customers a solution to a comfortable and enjoyable lifestyle even in the most unfortunate cases.

All mattresses are tested and certified by the Technical University of Berlin.





Bath and nursing appliances

Trautwein are experts in bath and nursing appliances since 1928. They have perfected their products in constant cooperation with nursing homes, hospitals and therapists.

Their in-house development department provides for the best available technology and security in our products. Innovations and many patents as well as a large number of clients stand for their high quality level.

The multifunctional product range is equipped for the bathing of people in need of care – and to spoil them. The special nursing tubs offer an extensive well-feeling experience.



THE HIGHEST QUALITY

The drive technology of the nursing tubs is tenthousandfold well proven. A very decisive point for Trautwein products is the high-grade material for a long-life cycle.

Trautwein uses high-quality fittings from Hans Grohe and produces all tub bodies with sanitary acrylic that is proven in many cases. The material is easy to clean and keeps its brilliant shine even after decades due to its hygienic level surface.

Free colour choice for tub body and cladding according to the Trautwein colour chart is offered at no extra charge.

TUB MATERIAL WITH ANTIBACTERIAL-TECHNOLOGY

Due to the antibacterial technology, this special acrylic tub material offers a durable and reliable protection against the growing of harmful bacteria like e.g. salmonella, E. coli, listeria, legionella, staphylococcus and streptococcus.



ThevoVital - The Mobility Bed

THEVOVITAL - DEMENTIA CARE

Breaking the vicious circle: Finally sleeping through the night despite dementia.

67% of patients with dementia suffer from sleep disturbances affecting both patient and caregivers. Sleep aid medications often become ineffective, and are not a long-term solution.



COMPONENTS

- + Mattress firm density foam provides a secure feeling when lying down
- + Wing suspensions adapt individually to the body zones
- + MiS Micro-Stimulation® improves body perception
- + Dementia suspension firm sensation of lying - good orientation is optimally tailored to the sleep and lying needs of persons with dementia
- + Improved sleep with the dementia mattress
- + Less night time wandering
- + Less napping during the day
- + Less disorientation
- + Decreased need for sleep aid medications
- + Peace of mind for caregivers
- + Documented experience reports
- + Two out of three patients with dementia showed a great improvement in their sleep patterns



HELPING AIDS

Bath Safety Quick-released shower chair without Back

- * Textured non-slip seat with drain holes
- * Light aluminum anodized frame
- * Durable and corrosion proof
- * 5 level height adjustment



Bath Safety Raised Toilet Seat 5"

- * Designed for individuals who have difficulty sitting down or standing up from toilet
- * Made of Aluminium tube
- * Fits most toilets



Walking Aid One Button Folding Walker

- * Aluminum anodized frame
- * 4 Legs have 8 level height adjustment



Walking Aid Delta Steel Rollator

- * Powder coated steel frame
- * Foldable design for travel
- * With basket, meal plate and bag



Push Button Drop Arm Commode with 4 Swivel Casters

- * Steel tubing welded construction
- * Bucket
- * Two lockable castors





HELPING AIDS

Aluminum Rollator

- * Aluminum liquid coated frame, with 1.6mm pipe thickness
- * Removable two-position basket with handle to be put in front of the rollator or under the seat
- * Easy squeeze and release hand brakes come with locking mechanism



Lightweight Wheelchair

- * Fixed footrests
- * Fixed full length armrests
- * Solid 200 x 25 front castors
- * Pneumatic 24" rear wheels



Standard Wheelchair

- * Fixed footrests
- * Fixed full length armrests
- * Solid 200 x 25 front castors
- * Solid 24" grey rear wheels



Commode

- * Height adjustable
- * Foldable
- * Weight capacity 90 kg
- * Safety frame over toilet



Locker with overbed table

- * Height adjustable tilting tray
- * Locker with door & drawer
- * 4 castors
- * Foldaway tray






PrimaCare

HOMECARE & CLINIC FURNITURE

Every person has a right to be treated with respect and dignity. Especially people in need of care. We think, develop and provide with this in mind.

Do not hesitate to contact us.

Hotline:

021 801 4675
071 206 6765

PrimaCare Homecare & Furniture (Pty) Limited
1 Printers Way, Unit 6, Montague Gardens
7447, Cape Town, South Africa